ExampleDocument

ISSUING WORK		R FROZEN CE	LLS AS STA				
		complete. Do not		-	-		
Patient Name				Date			
Patient Hospital No				DOB			
ID Number/s				Batch/N	umber		
Location R	ack	Slot		Date/s o	of Collection		
					Performed by	Checked	
Time of pickup confirm	ned			7			
Dry shipper received a	and cheo	ck list completed	d				
Check release authori	sed						
Date returned entered	in Syne	rgy					
Cassettes placed in a	shipper	to precool					
Rack set at 0 (if all ce	ls used)						
Tray prepared with dry	/ ice for	use checking th	ne bags				
Cryobag removed from Patient details, harves insert and product lab	t numbe	r and bag num		on			
Integrity of bag checke	ed						
Place in pre-cooled ca Cassette placed in po at least one absorber Polybag sealed and for cassette	ybags, (it sheet.	Cryobag ports	-	•			
Polybag placed in out	er bag,						
Outer bag sealed and the cassette	bag fold	led so that it wr	aps snugly a	round			
Bags packed individua to keep cool until pack	•	•		hipper			

ExampleDocument

Aliquots placed in inner leak-proof polybag with absorbent sheet, air removed and sealed.	
Aliquot bag placed in padded bag, air removed and sealed	
Check logger active and flashing green	
Rack removed and lid replaced	
Cassettes placed in rack Cryobag ports upwards	
Vials placed in rack	
Rack returned to shipper	
Time loaded into shipper Lid of inner shipper closed (holes must line up with holes in Dewar neck	
Serialised tag applied Tag number	
Outer shipper lid replaced	
Time placed in shipper	
Contents label attached to outer shipper	
Shipper certification form packed in outer pouch of shipper	
Latches sealed Copies of weigh bill (4) signed and placed in outer pouch(courier brings these).(Copy one for records)	
Handed to courier-courier signs transport form	
Time handed over	
Date and times recorded on Web Site and checked	