

SOP for preparation and maintenance of tumour-infiltrating lymphocyte, cancer associated fibroblast and autologous tumour cultures (ovarian cancer)

Lead Organisation: iMATCH

Author: The University of Manchester

Version Number: 1.0

Finalisation Date: 27/06/2025

End user rights:

This document is shared with permission for re-use to distribute, remix, adapt, and build upon the material in any medium or format for non-commercial purposes only, so long as the attributions listed below are given.

Attributions: University of Manchester, iMATCH collaboration, Manchester NHS Foundation Trust, Salford Royal NHS Trust, The Christie NHS Foundation Trust. Professor Richard Edmondson, Dominique Jones, Marcus Price.

This document is made available under a Creative Commons Attribution-NonCommercial 4.0 International License as described here: <https://creativecommons.org/licenses/by-nc/4.0/>

The information, materials and any opinions contained in this document are provided for general information and educational purposes only, are not intended to constitute legal, medical or other professional advice and should not be relied on or treated as a substitute for specific advice relevant to particular circumstances. Although we make all reasonable efforts to ensure the information is up to date, we make no representations, warranties or guarantees in that regard. In no event shall the creator(s) be liable for any direct, indirect, special, consequential or other claims, losses or damages that are related to the use or reliance whatsoever in the content of the document or any part thereof, except to the extent that such liability cannot be excluded by law. We do not seek to exclude or limit in any way our liability to the user for personal injury or death caused as a result of our negligence or seek to exclude or limit our liability for fraud or fraudulent misrepresentation by us.

We reserve the right to make changes and improvements to any information contained within this document, at any time and without notice. Where this document contains hyperlinks to other websites operated by parties not connected to us, such hyperlinks are provided for your reference only. We do not control such websites and are not responsible for their contents.

The inclusion of hyperlinks from this document or the website to such websites does not imply any endorsement of the material on such websites or any association with their operators. We accept no responsibility of any nature whatsoever for linked web sites or any information contained in them.

Funded by



Coordinated by



Standard Operating Procedure

(Version 1.1)

Preparation and maintenance of tumour-infiltrating lymphocyte, cancer associated fibroblast and autologous tumour cultures (ovarian cancer)

Revision History

Number	Date	Reason for Change
1.0		
1.1	12/09/2019	
1.2	06/12/2021	

Prepared By:

Checked By:

Authorised by:

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeonc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

**Gynaecological Oncology Group, Division of Cancer
 Sciences**

Signed: 	Signed:	Signed: Edmondson
Print: Caitlin Waddell	Print:	Print:
Date: 18.02.2019	Date:	Date: 21/2/22
Date of next review: 27.6.25	R. Edmondson No changes to content	Signed: R. Edmondson

*Uncontrolled document if photocopied and printed
 (Latest version available on Gynaeonc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

1. Introduction

Molecular and cell biology techniques are important tools for assessing the differences between normal and diseased cells and tissue. Such techniques have also been employed to detect markers that may predict some response to treatment and/or have prognostic significance. These translational research studies can greatly influence the success or failure of novel therapies. However, these investigations can only yield valid results if rigorous procedures are followed to ensure the collection and storage of high quality tissue samples.

2. Purpose

The purpose of this document is to provide clear guidance on the procedures to be followed during the processing of ovarian tumour tissue to establish a) a tumour infiltrating lymphocyte culture; (b) an autologous tumour cell culture; (c) a cryopreserved aliquot of tumour culture digest; (d) subsequent cryopreservation of TIL and autologous tumour cultures for future assays.

Departmental standard operating procedures will already exist within the institute/trust. These must be read and followed in conjunction with this document.

3. Location

This document pertains to processing of human ovarian tumour tissue on the 5th floor laboratories of St. Marys Hospital, Manchester.

4. Objective

The ultimate aims of this document are to ensure that TIL and tumour cultures are established reproducibly, and that appropriate records and recordings are taken to ensure that as much information as possible is available to allow analysis of data. Additionally, this document aims to ensure that cultures and disaggregates are cryopreserved to ensure that assays and/or rapid expansion protocol may be carried out on all patient's samples if needed.

5. Scope

This document refers to the collection of tissue for the MFT Biobank for use in the research carried out into ovarian cancer under Professor Richard Edmondson. The ultimate aims of this document are to ensure that tissue and fluid specimens, which are collected under the standard operating procedures established previously, are processed to produce cell cultures under the 'ovarian TIL' spectrum of research being carried out in the Edmondson group. Tissue collected on behalf of the lab will also be used to establish basic tumour culture (which will be cryopreserved and used for a raft of established characterisation assays), cryopreserved for RNA analysis, formalin fixed and paraffin embedded for immunohistochemistry, and some tissue gifted to collaborators under MTAs. Therefore, this SOP only concerns a small

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

portion of the work to be carried out on tissue procured through MFT Biobank, which should be considered in the process of carrying out this work.

6. Responsibilities

6.1. The collection of tissues and fluids containing malignant cancerous cells from patients suffering from ovarian cancer is the responsibility of the clinicians operating on or caring for the patient.

6.2. Following sampling of the tissues, the clinical team is responsible for contacting members of the Edmondson lab to ensure that samples are delivered and processed within 30 minutes of their sampling in theatre.

6.3. Clinical fellows and researchers are under the direct supervision of the institution / trust line management.

6.4 Clinical fellows and researchers, in collaboration with the clinical care team for the following:

6.5.1. Patient identification and consent, specimen labelling, early storage, and transportation of the samples

6.5.2. Ensuring that sampling and storage occurs within a timely manner
Once tissue is *ex vivo*.

6.5.3. Accurately completing any accompanying documentation

6.5.4. Investigation or being aware of back up procedures in the eventuality that a problem arises during sampling or transportation (e.g. blood spillage)

6.5.5 Liaising with other trust and Edmondson lab staff.

6.5 It is the responsibility of the staff or student processing procured tissue to ensure: (a) that tissue is processed in a timely manner; (b) that appropriate records are kept of tissue storage within the lab/storage facilities, and experimental procedures carried out on tissue and cell lines; (c) that

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

experiments are carried out safely and reproducibly, in accordance with guideline documents listed below; (d) that the staff/student member is familiar with the correct methods of disposal of tissue or downstream cultures/experimental products.

7. Related Documents

- 7.1. Guidelines for Good Clinical Practice (GCP)
- 7.2. Human Tissues Act 2004
- 7.3. COSHH regulations 1988
- 7.4. Local risk assessments / safe operating procedures
- 7.5. Local Standard operating procedures
- 7.6. National Blood service (www.blood.co.uk)
- 7.7. MFT Biobank list of SOPs and Guidance documents

Appendices

Appendix A: List of SOPs and Guidance documents

Standard operating procedures, Risk Assessments and related COSHH

- 'Disaggregation of solid tumour samples using gentleMACS dissociator and MACSmix tube rotator' (TIL1)
- 'Preparation of TIL Media and initiation and maintenance of TIL and autologous tumour cultures' (TIL2)
- 'Safe operation of a laboratory centrifuge'
- 'Use of a microbiological safety cabinet'
- 'Disposal of biological agents, human blood, body fluids and tissue'
- 'Primary cell culture from ascitic fluid and solid tumour tissue from patients undergoing treatment for ovarian cancer'

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

- 'Permitted laboratory disinfectants'
- 'How to deal with needle-stick injuries and contamination accidents'
- COSHH: Beta-mercaptoethanol (Gynae Onc drive)
- COSHH: Tumour Dissociation Kit, human: enzymes H, R and A (Gynae Onc drive)
- Further COSHH forms available for reagents used in this protocol are stored on the shared H&S drive.

Standard Operating Procedure

Tumour Disaggregation Protocol

Preparation:

- 1) Ensure that MACSmix rotator has been sufficiently charged to carry out incubation steps.
- 2) Warm 100mL RPMI-1640 (no supplements) and RPMI-1640 (+L-glutamine/+HEPES) for at least 2h before needed. Thaw any frozen aliquots or pen-strep stock. Avoid multiple freeze/thaws.
- 3) Prepare TIL media as below.

[C]	Reagent	Volume of stock (mL)	Stock [C]
	RPMI-1640 (supplements highlighted)	445	N/A
2mM	L-glutamine		2mM
25mM	HEPES		25mM

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeonc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

10%	HI FBS		50	N/A
100IU/mL	Penicillin	Combined reagent	5	10,000IU/mL
100µg/mL	Streptomycin			10,000µg/mL
50µM	β-mercaptoethanol		1µL/mL	50mM

4) Assemble instruments/consumables for tumour disaggregation:

- (a) Sterile scalpel and disposable forceps
- (b) 200mm petri dish
- (c) Balance
- (d) Miltenyi C-tubes
- (e) Tube racks (to fit 50mL Falcon/C-tube)
- (f) RPMI-1640 (no supplements)
- (g) Human tumour dissociation kit (Miltenyi Biotec)
- (h) P100/P200/P1000 and 20/200/1000µL tips

Protocol:

- 1) Weigh sample.
- 2) For up to 1g of tissue, add 4.7mL RPMI-1640 (no supplements) to the C-tube. Add 200µL enzyme H, 100µL enzyme R and 25µL enzyme A. For between 1 and 2g, double volumes of media and enzyme. Carefully secure the lid of the C-tube until it clicks.
- 3) Remove tumour sample from TSS. If sample is particularly bloody, it can be rinsed in fresh PBS. Remove any fat and weigh. Cut into 2-4mm³ fragments using a sterile scalpel. Rinse fragments again with PBS if needed. Place fragments into a C-tube (up to 2g/tube – run multiple tubes if necessary).
- 4) Follow dissociation program below using GentleMACS dissociator (GMD) and MACSmix (MmR) rotator in incubator:

Step No.	Equip.	Action	Time (mins:secs)	Temp (°C)
1	GMD	Run <i>h_tumour_01</i> program	0:36	RT

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

2	MmR	Rotate/incubate (20 rpm, cont. rot.)	30:00	37
3	GMD	Run <i>h_tumour_02</i> program	0:36	RT
4	MmR	Rotate/incubate (20 rpm, cont. rot.)	30:00	37
5	GMD	Run <i>h_tumour_03</i> program	0.36	RT

- 5) Use a 100µm cell strainer to filter digest into a 50mL Falcon. Wash filter with 20-40mL RPMI-1640. If digest is particularly fatty or fibrotic, the digest may need to be split between tubes/filters.
- 6) Centrifuge at 400xg for 7 minutes.
- 7) If processing omentum, gently aspirate fat layer from the surface of the supernatant with a 5-10mL stripette and discard before aspirating the remaining supernatant.
- 8) Aspirate supernatant and resuspend pellet in T-cell media (scale volume depending on pellet size). If digest is very contaminated with red blood cells, can use RBC lysis buffer (Miltenyi) before use (mandatory if being used for flow).
- 9) Count live cells using Trypan blue.

Initiation of Tumour/CAF Cultures

- 1) Transfer desired volume of tumour digest to a fresh tube and centrifuge at 400xg for 5 minutes.
- 2) Aspirate supernatant and resuspend pellet in RPMI 1640 with 10% FBS, 1% penicillin-streptomycin and 0.1% amphotericin-B.
- 3) Transfer to a 50µg/cm² collagen-I coated vessel (approximate plating densities <math><0.5 \times 10^6</math> cells into T25, - 4) After 24 hours, assess adherence and confluency.
- 5) After 48-72 hours, dead/unattached cells can be washed off with DPBS and media replaced.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

- 6) When vessel reaches approx. 70% confluency, assess morphology of cells (look for colonies of cobblestone-like cells).
- 7) Wash flask with DPBS and add 2-4mL of 0.05% trypsin-EDTA.
- 8) Assess detachment of cells under light microscope throughout, gently swirling the trypsin over the adherent surface every minute. 'Fibroblasts' should round up quickly (<2 mins) while cobblestone-like colonies will remain attached for much longer.
- 9) If detachment after rounding up is slow, the flask can be placed in the incubator for a minute before checking again.
- 10) When the majority of cells surrounding cobblestone colonies have detached, aspirate trypsin gently and add to 10mL of DMEM with 10% FBS, 1% penicillin-streptomycin and 0.1% amphotericin.
- 11) Wash the adherent surface of the flask very gently with DPBS (add DPBS, swirl flask and aspirate rather than wash the surface directly with a stripette) and add to the DMEM 'fibroblast' tube.
- 12) Add complete RPMI back to the trypsinised flask as soon as possible and return to hypoxic incubator.
- 13) Centrifuge the fibroblasts at 400xg for 5 mins.
- 14) Resuspend pellet in complete DMEM and count.
- 15) Transfer to 0.1% gelatin coated flask (follow approximate density guidelines in point 3) and incubate in 5% CO₂ incubator.
- 16) Wash and change media of both lines every 2-3 days, sub-culturing at 70% confluency.

Initiation of TIL Cultures

Preparation:

- 1) Using cell count from Trypan blue step: calculate volume of media needed to achieve a final cell density of 0.5x10⁶ cells/mL (2mL) in 24 well plates (TC treated).
e.g. Digest has yielded x cells/mL, need 0.5x10⁶ cells/mL.

$$\text{Volume of media} = \frac{\text{digest cell count (cells/mL)} \times \text{volume}}{0.5 \times 10^6 \text{ cells/mL}}$$

As 2mL/well is required for the 24 well plate, divide the volume by two to work out how many wells and plates can be filled.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeonc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

- 2) Based on volume calculation, determine the volume of IL-2 (once vial reconstituted, stored in -80°C at a stock concentration of 1MIU/mL) needed to achieve 3000IU/mL.
- 3) Finally, based on 0.5×10^6 cells/mL, calculate the volume of Dynabeads (Human T-Activator CD3/CD28) needed to achieve a 1:1 ratio of Dynabeads:cells. Dynabead concentration = 4×10^7 beads/mL.

$$\text{Volume of Dynabeads} = \frac{0.5 \times 10^6 \text{ cells/mL}}{4 \times 10^7 \text{ beads/mL}} = 12.5 \mu\text{L Dynabeads/mL}$$

- 4) Prepare Dynabeads for use:
 - (a) Resuspend Dynabeads within the vial by vortexing for >30 seconds, and then transfer the required final volume to a 15mL Falcon.
 - (b) Add 1mL of PBS, and vortex for 5-10 seconds.
 - (c) Place tube on a magnet until all beads are captured, and discard the supernatant.
 - (d) Resuspend beads in the same volume of medium as original Dynabead volume.

Protocol:

- 1) Resuspend cell pellet in calculated T-cell media volume to final cell density of 0.5×10^6 cells/mL.
- 2) Add 12.5 $\mu\text{L/mL}$ Dynabeads to cell suspension. Next, add 3000IU/mL IL-2. Invert gently and thoroughly to mix.
- 3) Split suspension between wells of 24 well plates, using 2mL/well.
- 4) Label plate with MOC number, tissue type, plate number, conditions, date, investigator.
- 5) Incubate at 37°C, 5% CO₂ and 95% humidity. Follow protocol below for growth phase.

Day	Action	Day	Action
0	Seed plates.	13	Count, reseed.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Gynaecological Oncology Group, Division of Cancer Sciences

2	Half media change (3000IU/mL IL-2)	15	Count, reseed.
4	Half media change (3000IU/mL IL-2)	17	Count, reseed.
6	Half media change (3000IU/mL IL-2)	19	Assays, freeze down, rapid expansion.
7	Remove Dynabeads, count, reseed.		
9	Count, reseed.		
11	Count, reseed.		

Half media change (day 2-6):

Note: TILs are in suspension, so will settle at the bottom of the plate but will not be stuck down. Plates must be moved very carefully from the incubator, and media aspirated only from the very top of the media layer. Do not attempt to aspirate media if plates have been disturbed.

- 1) Make up volume of T-cell media with 6000IU/mL IL-2.
- 2) View wells on microscope, check for signs of infection (cloudy media, yellowing etc.).
- 3) Carefully remove 1mL media from each well using a P1000 and barrier tip. Discard media into Virkon.
- 4) Slowly add supplemented T-cell media to each well and return to incubator.

Dynabead removal (day 7):

- 1) Thoroughly resuspend wells and add to Falcon tube. Rinse wells with T-cell media.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeonc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

- 2) Remove Dynabeads using Invitrogen magnet. Allow 5 minutes to ensure that all beads are captured.
- 3) Remove supernatant to fresh tube.
- 4) Wash beads by resuspending in fresh media and capture beads again. Add supernatant to first harvest of cells.
- 5) Spin down at 400xg for 5 mins, count using Trypan blue and reseed in 24 well plates or flask at 1×10^6 cells/mL in T-cell media supplemented with 1000IU IL-2.

$$\text{Volume of media} = \frac{\text{cell density (cells/mL)} \times \text{volume}}{1 \times 10^6 \text{ cells/mL}}$$

Count and reseeding (day 9-17): Repeat steps 1 and 5, reseeding in T-cell media with 1000IU/mL IL-2 to a final density of 1×10^6 cells/mL.

Rapid Expansion Protocol

Preparation:

- 1) Calculate volume of TIL suspension containing 1×10^5 TIL.
- 2) Calculate volume of irradiated feeder cell suspension (30Gy) containing 2×10^7 cells. Need a 1:200 TIL:feeder cell ratio.
- 3) Prepare stock solution of anti-CD3 Ab (OKT3) if necessary. Calculate volume of Ab needed for final concentration of 30ng/mL in 20mL media.

Protocol:

- 1) Seed a T75 flask with correct volumes of TIL suspension and feeder cells, suspended in 20mL T-cell media supplemented with 30ng/mL OKT3 and 3000IU/mL IL-2.
- 2) Place in incubator. Leave undisturbed for 3-5 days, until cells start to colonise.
- 3) Then replace media every 3-5 days, until flask reaches 80% confluency, and subculture.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Tumour Culture Initiation Protocol

Preparation:

- 1) Warm 500mL RPMI-1640 (+L-glutamine, +HEPES) and 1X PBS, along with any frozen aliquots for at least 45 minutes before they are to be used in culture.
- 2) Prepare tumour cell 'complete' media as in table below:

[C]	Reagent		Volume of stock (mL)	Stock [C]
	RPMI-1640 (supplements highlighted)		395	N/A
2mM	L-glutamine			2mM
25mM	HEPES			25mM
20%	HI FBS		100	N/A
100IU/mL	Penicillin	Combined reagent	5	10,000IU/mL
100µg/mL	Streptomycin			10,000µg/mL

Protocol:

- 1) Transfer digested tumour from C-tube to 50mL Falcon and spin in centrifuge at 400xg for 5 mins.
- 2) Discard the supernatant and wash the pellet with 1X PBS 3 times.
- 3) After the third wash, resuspend the pellet in 'complete' media. Count the cells, and seed in a T75 flask to a volume of 20mL.
- 4) Label the flask with MOC number, tissue type, date, investigator, flask number.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)

Gynaecological Oncology Group, Division of Cancer Sciences

- 5) Place in top incubator and leave undisturbed for 3-5 days, until cells are approximately 80% confluent. If media becomes discoloured, replace half of the media.
- 6) Upon confluency, passage cells as below.

Subculture/passage of tumour cells:

- 1) Remove old media. Wash cells once with 1X PBS, then detach using enzyme-free cell dissociation buffer (in PBS).
- 2) Incubate for 2-5 mins at RT. Aspirate into a 50mL Falcon and add an equal volume of complete medium to dilute enzyme.
- 3) Centrifuge cells at 400xg for 5 minutes.
- 4) Remove supernatant and resuspend in 30mL of fresh media.
- 5) Transfer cells into 3x T75, and suspend in 10-15mL fresh media.
- 6) Label flask with MOC number, tissue type, date, investigator, flask number.

Cryopreservation of TIL/ autologous tumour cultures and tumour disaggregate

- 1) Prepare freezing media using 90% HI-FBS and 10% DMSO.
- 2) Transfer relevant cells to a 50mL Falcon. Count cell density.
- 3) Centrifuge at 400xg for 5 mins to form a cell pellet, discard supernatant, and resuspend in appropriate volume of freezing media (see below).

Cell type	Cell density (cells/mL)
Tumour digest	>2x10 ⁷
Expanded TIL	>2x10 ⁷
Cell lines	>1x10 ⁶

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeconc shared drive)*

**Gynaecological Oncology Group, Division of Cancer
Sciences**

Tumour culture	>2x10 ⁷
----------------	--------------------

- 4) Prepare 1mL aliquots of cells in freezing media in 1.6mL cryovials, and label with MOC number, cell type, passage, cell density, date and initials.
- 5) Place in -80°C storage in a Mr. Frosty. Transfer to LN₂ dewar for long term storage.

*Uncontrolled document if photocopied and printed
(Latest version available on Gynaeonc shared drive)*

Preparation and maintenance of tumour-infiltrating lymphocyte and autologous tumour cultures (ovarian cancer)